



LexaMed 705 Front Street, Toledo, OH 43605 Phone 419-693-5307 Fax 419-691-0418

Growth-Chek™ Microbial Suspensions Certificate of Performance

Product: *Staphylococcus aureus*
Dose: < 100 CFUs / 0.1 mL
Passage Number: 4
Product Number: 6538-100
Lot Number: 031218
Expiration Date: June 12, 2018

Each vial contains <100 viable organisms per 0.5mL or 0.1mL dose. Each lot is derived from cultures traceable to a recognizable culture collection identified in the current *United States Pharmacopeia* (USP) and is certified for purity, and population and is not more than five passages from the original culture.

Certified by: Yamara Aller **Date:** 03/19/18
QA Approval: David R. Hakman **Date:** 3-19-18