



LexaMed 705 Front Street, Toledo, OH 43605 Phone 419-693-5307 Fax 419-691-0418

Growth-Chek™ Microbial Suspensions Certificate of Performance

Product: *Kocuria rhizophila*
Dose: < 100 CFUs / 0.1 mL
Passage Number: 5
Product Number: 9341-100
Lot Number: 102317
Expiration Date: December 23, 2017

Each vial contains <100 viable organisms per 0.5mL or 0.1mL dose. Each lot is derived from cultures traceable to a recognizable culture collection identified in the current *United States Pharmacopeia* (USP) and is certified for purity, and population and is not more than five passages from the original culture.

Certified by: *Samara Allen* **Date:** 10/30/17

QA Approval: *BA* **Date:** ~~11-3~~ 11-1-17
11-1-17 BA wrong date