



LexaMed 705 Front Street, Toledo, OH 43605 Phone 419-693-5307 Fax 419-691-0418

Growth-Chek™ Microbial Suspensions Certificate of Performance

Product: *Kocuria rhizophila*
Dose: < 100 CFUs / 0.5 mL
Passage Number: 5
Product Number: 9341-20
Lot Number: 051418
Expiration Date: July 14, 2018

Each vial contains <100 viable organisms per 0.5mL or 0.1mL dose. Each lot is derived from cultures traceable to a recognizable culture collection identified in the current *United States Pharmacopeia* (USP) and is certified for purity, and population and is not more than five passages from the original culture.

Certified by: *Jamara Allen* Date: *05/23/18*
QA Approval: *Brand L. Heckman* Date: *5-23-18*