



*LexaMed* 705 Front Street, Toledo, OH 43605 Phone 419-693-5307 Fax 419-691-0418

---

## Growth-Chek™ Microbial Suspensions Certificate of Performance

Product: *Kocuria rhizophila*  
Dose: < 100 CFUs / 0.5 mL  
Passage Number: 5  
Product Number: 9341-20  
Lot Number: 061818  
Expiration Date: August 18, 2018

Each vial contains <100 viable organisms per 0.5mL or 0.1mL dose. Each lot is derived from cultures traceable to a recognizable culture collection identified in the current *United States Pharmacopeia* (USP) and is certified for purity, and population and is not more than five passages from the original culture.

Certified by: Jamara Allen Date: 06/27/18  
QA Approval: Brandi K. Heck Date: 6-27-18